OZONFIX^{oliv}

SKŁAD:

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OZONFIX oliv

COMPOSITION:

Fatty acids: oleic acid, linolenic acid, palmitic acid, stearic acid, α -linolenic acid and oligomeric ozonides of these acids (polyperoxides), oxygen ions not less than 4,9 g/l.

PROPERTIES

OZONFIX oliv has very similar biocidal properties to pure ozone.

This product meets the standards of PN-EN 1040:2006: Chemical disinfectants and antiseptics.

OZONFIX oliv is STERILE.

It is a natural sunscreen, SPF=7,5

This product is perfectly absorbed into the skin, transporting oxygen ions into the dermis.

It oxidizes histamine and free hemoglobin, oxygenates tissues at the site of application. Active oxygen is directly involved in cellular reactions.

It accelerates metabolism at the site of application, it may be used for virus inactivation and for killing pathogenic microorganisms such as bacteria and fungi. OZONFIX oliv soothes skin irritation, protects damaged skin (mucosa) from bacterial and fungal complications.

It eliminates pain and itching.

The released oxygen ions stimulate the production of cytokines, which pro-

motes wound granulation. This product makes the epidermis more flexible.

APPLICATION

This product can ONLY be used for external use on the skin, nails and mucous membranes. There are no age restrictions.

OZONFIX oliv reduces the unpleasant effects of, for example, bee, mosquito and blackfly bites. It has anti-allergic and anti-inflammatory effects.

It can be applied in the case of diaper dermatitis, sunburn, chemical burn, irradiation with the cobalt bomb, frostbite. Oxygenated tissues revitalise faster. In the case of acne, they relatively reduce the number of pimples due to faster skin revitalization and protection of the pimples against bacterial complications (when used regularly).

For those with atopic skin, psoriasis (when used as a supportive product), eczema, diabetic foot syndrome or bedsores.

OZONFIX oliv perfectly moisturizes the skin and contained oxygen ions protect the skin from the external microenvironment. The skin is oxygenated. In the case of bedsores and diabetic foot, OZONFIX oliv should be used as an oxygenator (also preventively).

This product can be applied to fresh post-operative wounds – they revitalise faster.

Thanks to instant killing of viruses, such as herpesviruses, this product prevents from the formation of lesions on the lips. If they already are, OZONFIX oliv accelerates lip revitalization. This product can work for chickenpox in a similar way; it additionally eliminates itching. Apply OZONFIX oliv in the case of oral mucositis, oral candidiasis, fungal infection (also on genital area).

Because of its antibacterial properties, this product can also be applied in the case of, for example, otitis, complications of wounds.

OZONFIX oliv can be used for cosmetic purposes, firming the skin, skin oxygenation (oxygen infusion), fighting cellulite. This product accelerates the absorption of subcutaneous hematoma and petechiae (apply OZONFIX oliv immediately after the incident, as it does not work on existing bruises).

OZONFIX oliv is sterile. Therefore it can be applied to open wounds.

This product can be used as a disinfectant after the removal of hybrid manicure from nails, cuticle cutting, depilation and after other invasive body treatments (peeling, mesotherapy, Plasma Complex, etc.) – OZONFIX oliv relieves pain and accelerates revitalisation, protects against bacterial and fungal complications.

OZONFIX oliv is a safe product for everyday use because it DOES NOT MAKE MICROORGANISMS RESISTANT TO ITS EFFECTS.

APPLICATION METHODS

For disinfection use, OZONFIX oliv needs to be applied in a sterile environment using the OZONFIX oliv pipette

dropper (included in the package) on the wound and its surrounding area. A hint for the product to work better: the wound should be clean and free of necrotic tissue. Apply this product 2 to 3 times daily for at least 5 days. In the case of shallow wounds, there is no need to apply a dressing – just rub OZONFIX oliv lightly. For deep wounds, apply wet dressings by soaking a sterile gauze in OZONFIX oliv, leaving it on for longer time. Change dressings at least once a day. OZONFIX oliv can be used for cleansing wounds.

Cosmetically OZONFIX oliv works better when applied to dry skin.

SIDE EFFECTS, CONTRAINDICATIONS

This product has an antihistamine effect, even through the skin. Redness after the application of OZONFIX oliv may indicate a strong reaction with oxidised histamine. The redness should subside after several minutes without any feeling of discomfort.

Ingestion of OZONFIX oliv results in depletion of the gastrointestinal tract and diarrhoea. Probiotics should be used.

STORAGE

Store OZONFIX oliv away from light and water. Keep it at room temperature in a hermetically sealed package. Shelf life of the product – 12 months after opening of the package.

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